

SWIFT

WWW.ROTAID.COM



 **ROTAID™**
AED CABINETS

ROTAID SWIFT

PRODUCT SHEET

The Rotaïd Swift is an innovative AED wall cabinet that ensures a recognizable and sustainable storage of your AED device. The characteristic shape and green color make it easy to be found and recognized and comes with a high quality appearance. A quick twist to the left opens the cabinet in an instant.

The large window, with a prominent AED logo, offers a full view on the defibrillator and its indicator.



An optional **audible alarm** of 90 dB makes sure the AED is well protected.

An **adjustable stand** allows placement of **every type of AED**.

The **solid steel back plate** enables installation on any surface and offers a high quality appearance.

The Swift is opened by an **intuitive twist** that removes the whole cover and provides the responder with full access to the AED.

Dimensions

Length **43,5 cm** Width **36,6 cm** Depth **21 cm** Weight **4,5 kg**

Length **17,1"** Width **14,4"** Depth **8,3"** Weight **9,9 lbs**